**ANCHIPLAST s.r.l.**

Via Briana, 4 I-46023 BONDENO DI GONZAGA (MN) ITALY

(REF:6400,6401,6404,6405,6407,6409,6410,6415,6416,6417,6419,6423,6427,6434,6435,6439,6444,6445,6446,6448,6449,6452,6453,6455,6456,6458,6460,6461,6462)

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#### **SET FOR DENTAL IRRIGATION**

Disposable monouse and single-patient device. Do not reuse and dispose of, after its use, according to the regulation in force. Its reuse implies risk of infections to the patient and/or user and can compromise the correct functionality of the device. Check that the device is compatible with the equipment with which it is to be used. Check as stated in the name of the device.

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#### **PHTHALATE FREE PRODUCT INSTRUCTION FOR USE**

1. Disinfect bottle cap, insert the piercing and overturn the bottle. 2. Close the flow regulator, press the tank until it's half full. 3. Insert the silicon attachment around the pump and secure the plastic locks in the slot (please refer to the unit use manual). 4. Connect the tube to the contrangle. 5. Open the flow regulator. Start the pump and completely fill the tube. Let little solution come out.

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**STERILE, OPERATE UNDER ASEPTIC CONDITIONS, DO NOT USE IF PACKAGE IS DAMAGED, WE ADVISE NOT OVERCOME THE MAXIMUM LIMIT OF 0,5 bar.**

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1. Connect the end of the irrigation tubing onto the handpiece irrigation nozzle. 2. Use the 5 clips provided in order to clip together the irrigation tubing and handpiece cable. 3. Insert the thicker tubing segment of the sterile line into the pump. 4. Remove the protective cap from the spike and connect it to the irrigation bag. 5. Open the air inlet on the tube before operating. Release the irrigation tube clamp, if this is closed. 6. Warning: Single use only. Disposed of the used irrigation tubing in accordance with national or local regulations or guidelines, as appropriate.

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**Check that the device is compatible with the equipment with which it is to be used. Check as stated in the name of the device.**

**WARNING: There are no specific contraindications related to the use of the device. Contraindications related to dental irrigation exists, therefore a careful monitoring of the patients during treatment and the evaluation of pertinent medical literature is recommended. Treatment shall be interrupted immediately in case of complications affecting patient's conditions. Stop using the device in case of allergic reactions to any component.**

**WARNING: Notice to the user and/or patient: report any serious incident occurred in relation to the device to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.**

**WARNING: Use of the device is reserved to authorized medical staff only.**

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#### **SET FOR DENTAL IRRIGATION**

Disposable and single-patient device. Do not reuse and dispose of, after its use, according to the regulation in force. Its reuse implies risk of infections to the patient and/or user and can compromise the correct functionality of the device. Check that the device is compatible with the equipment with which it is to be used. Check as stated in the name of the device.

STERILE, OPERATE UNDER ASEPTIC CONDITIONS, DO NOT USE IF PACKAGE IS DAMAGED, WE ADVISE NOT OVERCOME THE MAXIMUM LIMIT OF 2 bar.

#### **PTHALATE FREE PRODUCT, NOT MADE WITH NATURAL RUBBER LATEX INSTRUCTION FOR USE**

1. Connect the end of the irrigation tubing onto the handpiece irrigation nozzle. 2. Use the 5 clips provided in order to clip together the irrigation tubing and handpiece cable. 3. Insert the thicker and 15 cm long tubing segment of the sterile line into the pump. 4. Remove the protective cap from the spike and connect it to the irrigation bag. 5. Open the air inlet on the tube before operating. Release the irrigation tube clamp, if this is closed. 6. Warning: Single-use only. Disposed of the used irrigation tubing in accordance with national or local regulations or guidelines, as appropriate.

**WARNING:** There are no specific contraindications related to the use of the device. Contraindications related to dental irrigation exists, therefore a careful monitoring of the patients during treatment and the evaluation of pertinent medical literature is recommended. Treatment shall be interrupted immediately in case of complications affecting patient's conditions. Stop using the device in case of allergic reactions to any component.

**WARNING:** Notice to the user and/or patient: report any serious incident occurred in relation to the device to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

**WARNING:** Use of the device is reserved to authorized medical staff only.

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